

Standard transformation using heat-shock method

🍏 Frozen competent cell preparation

- Overnight culture *E. coli* in SOB medium (eg. DH5 α) (5ml)
- Subculture to 50 ml SOB medium (1:100), shake at 37°C, 30 mins (O.D.=0.3)
- Transfer culture into centrifuge tube, put on ice for 10-15 mins
- Centrifuge 2000-3000rpm for 12-15 mins at 4°C
- Collect pellet, and resuspend in 16.6ml FSB, put on ice for 10-15 mins
- Centrifuge 2000-3000rpm, 12-15 mins at 4°C
- Resuspend in 4ml FSB
- + 140 μ l DMSO (i.e. to 3.5%), put on ice for 10 mins
- + 140 μ l DMSO (i.e. to 7%), put on ice for 10-20 mins
- Quick freezing by placing on dry-ice or liquid nitrogen, then store at -70°C

🍏 Transformation

- Take 100 μ l (or 200) competent cell, + 5-10 μ l ligased DNA
- Put on ice for 30 mins
- Heat shock at 42°C, 2 mins
- Put back on ice for 2 more mins
- Add 400 μ l SOC medium (or LB) per tube
- Incubate and shake at 37°C, 30-60 mins
- Spread the cells to plates (50 μ l, 200 μ l)

🍏 SOB medium (for 500 ml):

tryptone	2%	10g
yeast	0.5%	2.5g
NaCl	10mM	0.3g (5M, 1ml)
KCl	2.5mM	0.093g (1M, 1.25ml)
MgCl ₂	10mM	1M, 5ml
MgSO ₄	10mM	1M, 5ml

🍏 SOC medium: SOB + 20mM glucose