

RNA extraction method*

0.8g of tissue

Jer-Ming Hu Nov. 2000

Before you do, check the needed items (next page)

- 1. Add 8ml TLE, 800 μ l 10% SDS and 2.7ml of phenol to a 30ml Corex tube. Put on ice.**
- 2. Chill mortar and pestle with liquid nitrogen. Grind 0.8g tissue to fine powder and add to previous Corex tube with scratch rods.**
Place a plastic weigh boat floating on liquid nitrogen and sweep the powder into the boat.
- 3. Put in Fisher tissue homogenizer and mangle for ca. 2 mins.**
Put other samples on ice and also all following steps.
- 4. Add 2.7ml chloroform and homogenize for another 2 mins.**
- 5. Cover tube w/ parafilm and centrifuge at 12,000g, 4°C for 20mins.**
- 6. Transfer upper layer to clean tube and add 5ml of phenol/chloroform (1:1). Centrifuge at 12,000g, 4°C for 20mins.**
- 7. Repeat step 6.**
- 8. Transfer upper layer to clean tube and add 5ml of chloroform and centrifuge as above.**
- 9. Transfer upper layer to a very clean 15ml Corex tube and add 1/3 volume 8M LiCl. Leave at cold room overnight.**
Alternatively you can continue to next step immediately if you see lots of precipitation.
- 10. Centrifuge at 14,000g for 20mins, 4°C. Do a quick 70%EtOH wash.**
- 11. Dry the pellet and redissolve in 330 μ l of DEPC-H₂O, stay on ice; transfer to a clean 1.7ml eppendorf. Add 110 μ l of 8M LiCl. Vortex and sit on ice for >2hrs.**
If the pellet is very difficult to resuspend, add double amount of water and LiCl.
- 12. Centrifuge at 14,000g, 20mins, 4°C. Remove supernatant, wash with 70% EtOH briefly, dry and resuspend in 125 μ l DEPC-H₂O. Add 12.5 μ l 3M NaOAc and 355 μ l 95% EtOH. Vortex and sit for >2hrs at -20°C.**
- 13. Centrifuge at 14,000g, 20mins, 4°C. Remove supernatant, wash w/ 70% EtOH briefly, dry and resuspend in 50 μ l DEPC-H₂O. Take OD reading and/or run on formaldehyde gel w/ standards.**

*This method is based on the one described by Michael Frohlich 3/99, which was modified from Steve Jacobsen's and Neil Olszewski's methods.

The items you need for clean and prepared:

Chloroform washed 30ml (4/sample) and 15ml (1/sample) COREX tubes
Mortar and pestle
New or chloroform washed scratch rods
DEPC treated water
TLE (0.2M Tris, 0.1M LiCl, 5mM EDTA, pH8.2) (RNase-free)
10% SDS (RNase-free)
8M LiCl (RNase-free)
Phenol (~8ml/sample); chloroform (13ml/sample)
70% EtOH (RNase-free)
95% EtOH (RNase-free)
3M NaOAc (RNase-free)
Liquid Nitrogen