

**PRODUCT INFORMATION**

**Terminal Deoxynucleotidyl Transferase**

Pub. No. MAN0013724

Rev. Date 30 September 2016 (A.00)

#\_

Lot: \_ Expiry Date: \_

**Store at -20 °C**

Components	#EP0161	#EP0162
Terminal Deoxynucleotidyl Transferase, 20 U/μL	500 U	2500 U
5X Reaction Buffer	0.4 mL	2 × 1 mL



**DANGER**

Toxic if inhaled. May cause cancer. Toxic to aquatic life with long lasting effects. Avoid breathing dust/fume/gas/mist/vapours/spray. Use personal protective equipment as required. IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Dispose of contents/container in accordance with local/regional/national/international regulations.

Thermo Fisher Scientific Baltics UAB, V.A. Graiciuno 8, LT-02241 Vilnius, Lithuania  
tel. +370 700 55131 Note. For more information please see the Product Information.

[www.thermofisher.com](http://www.thermofisher.com)

**For Research Use Only.** Not for use in diagnostic procedures.

**Description**

Terminal Deoxynucleotidyl Transferase (TdT), a template-independent DNA polymerase, catalyzes the repetitive addition of deoxyribonucleotides to the 3'-OH of oligodeoxyribonucleotides and single-stranded, or double-stranded DNA (1). The TdT requires an oligonucleotide of at least three nucleotides to serve as a primer. With RNA as template TdT shows variable performance which strongly depends upon the tertiary structure of acceptor RNA 3'-end and the nature of nucleotide. Generally, it is lower than using DNA as a template.

**Applications**

- Production of synthetic homo- and heteropolymers (1).
- Homopolymeric tailing of linear duplex DNA with any type of 3'-OH terminus (2, 3), see protocol on back page.
- Oligodeoxyribonucleotide and DNA labeling (2, 4-8), see protocol on back page.
- 5'-RACE (Rapid Amplification of cDNA Ends) (9).
- *In situ* localization of apoptosis (10).

**Source**

*E.coli* cells carrying a cloned gene encoding calf thymus terminal deoxynucleotidyl transferase.

## Definition of Activity Unit

One unit of the enzyme catalyzes the incorporation of 1 nmol of deoxythymidylate into a polynucleotide fraction in 60 min at 37 °C.

## Storage Buffer

The enzyme is supplied in: 100 mM potassium acetate (pH 6.8), 2 mM 2-mercaptoethanol, 0.01% (v/v) Triton X-100 and 50% (v/v) glycerol.

## 5X Reaction Buffer

1 M potassium cacodylate, 0.125 M Tris, 0.05% (v/v) Triton X-100, 5 mM CoCl<sub>2</sub> (pH 7.2 at 25 °C).

## Inhibition and Inactivation

- Inhibitors: metal chelators, ammonium, chloride, iodide, phosphate ions.
- Inactivated by heating at 70 °C for 10 min or by addition of EDTA.

## Note

Due to the presence of CoCl<sub>2</sub> the TdT Reaction Buffer is incompatible with downstream applications. It is necessary to remove CoCl<sub>2</sub> from the reaction mixture by spin column or phenol/chloroform extraction and subsequent ethanol precipitation.

## CERTIFICATE OF ANALYSIS

### Endodeoxyribonuclease Assay

No detectable degradation was observed after incubation of supercoiled plasmid DNA with Terminal Deoxynucleotidyl Transferase.

### Ribonuclease Assay

No detectable degradation was observed after incubation of [3H]-RNA with Terminal Deoxynucleotidyl Transferase.

### Labeled Oligonucleotide (LO) Assay

No detectable degradation after incubation of single stranded and double stranded radiolabeled oligonucleotides with Terminal Deoxynucleotidyl Transferase.

Quality authorized by:



Jurgita Zilinskiene

*(continued on back page)*

## Protocol for tailing of DNA 3'-termini

1. Prepare the following reaction mixture:

<b>5X reaction buffer for Terminal Deoxynucleotidyl Transferase</b>	4 $\mu$ L
<b>DNA fragments</b>	1 pmol of 3'-ends
<b>dATP or dTTP or dGTP or dCTP</b>	130 pmol or 60 pmol
<b>Terminal Deoxynucleotidyl Transferase</b>	1.5 $\mu$ L (30 U)
<b>Water, nuclease-free (#R0581)</b>	to 20 $\mu$ L

2. Incubate the mixture at 37 °C for 15 min.

3. Stop the reaction by heating at 70 °C for 10 min or by the addition of 2  $\mu$ L 0.5 M EDTA (#R1021).

### Note

- Under the conditions described above, 100-130 dA or dT residues, or 20-30 dC or dG residues can be added per 3'-OH end of DNA.
- The efficiency of the reaction depends upon the type of 3'-OH termini of the DNA fragments. 3'-overhangs are tailed with higher efficiency than recessed or blunt ends.

## Protocol for DNA and oligonucleotide 3'-end labeling by tailing

1. Prepare the following reaction mixture:

<b>5X reaction buffer for TdT</b>	10 $\mu$ L
<b>Linear DNA</b>	10 pmol
<b>[<math>\alpha</math>-<sup>32</sup>P]-ddATP, ~10 TBq/mmol (3000 Ci/mmol)</b>	1.85 MBq (50 $\mu$ Ci)
<b>Terminal Deoxynucleotidyl Transferase</b>	2 $\mu$ L (40 U)
<b>Water, nuclease-free (#R0581)</b>	to 50 $\mu$ L

2. Incubate the mixture at 37 °C for 15 min.

3. Stop the reaction by heating at 70 °C for 10 min or by adding 5  $\mu$ L 0.5 M EDTA (#R1021).

### Note

The efficiency of the reaction depends upon the type of 3'-OH termini of the DNA fragments. 3'-protruding ends are labeled with higher efficiency than recessed or blunt ends.

## References

1. Bollum, F.J., Terminal deoxynucleotidyl transferase, The Enzymes, the third edition (Boyer, P.D., ed.), Academic Press, New York, vol.10, 145-171, 1974.
2. Deng, G.R., Wu, R., Terminal transferase: Use in the tailing of DNA and for *in vitro* mutagenesis, Meth. Enzymol., 100, 96-116, 1983.
3. Eschenfeldt, W.H., et al., Homopolymeric tailing, Meth. Enzymol., 152, 337-342, 1987.
4. Tu, C.-P.D., Cohen, S.N., 3'-end labeling of DNA with [ $\alpha$ - $^{32}$ P]-cordycepin-5'-triphosphate, Gene, 10, 177-183, 1980.
5. Vincent, C., et al., Synthesis of 8-(2,4-dinitrophenyl-2,6-aminohexyl)aminoadenosine-5'-triphosphate: Biological properties and potential uses, Nucleic Acids Res., 10, 6787-6796, 1982.
6. Kumar, A., et al., Nonradioactive labeling of synthetic oligonucleotide probes with terminal deoxynucleotidyl transferase, Anal. Biochem., 169, 376-382, 1988.
7. Gaastra, W., Klemm, P., Radiolabeling of DNA with terminal transferase, Methods in Molecular Biology, vol.2: Nucleic Acids (Walker, J.M., ed.), Humana, Clifton, NJ, 269-271, 1984.
8. Igloi, G.L., Schiefermayr, E., Enzymatic addition of fluorescein- or biotin-riboUTP to oligonucleotides results in primers suitable for DNA sequencing and PCR, BioTechniques, 15, 486-497, 1993.
9. Frohman, M.A., et al., Rapid Production of Full-Length cDNAs from Rare Transcripts: Amplification Using a Single Gene-Specific Oligonucleotide Primer. Proc. Natl. Acad. Sci. USA, 85, 8998-9002, 1988.
10. Gorczyca, W., et al., Detection of DNA strand breaks in individual apoptotic cells by the *in situ* terminal deoxynucleotidyl transferase and nick translation assays, Cancer Res., 53, 1945-1951, 1993.

### **LIMITED USE LABEL LICENSE: Internal Research and Development Use Only.**

The purchase of this product conveys to the buyer the limited, non-exclusive, non-transferable right (without the right to resell, repackage, or further sublicense) to use this product for internal research and development purposes. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of the product does not include or carry any right or license to use, develop, or otherwise exploit this product commercially and no rights are conveyed to the buyer to use the product or components of the product for purposes including but not limited to provision of services to a third party, generation of commercial databases or clinical diagnostics. This product is sold pursuant to authorization from Thermo Fisher Scientific and Thermo Fisher Scientific reserves all other rights. For information on purchasing a license for uses other than internal research and development purposes, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies Inc., 5791 Van Allen Way, Carlsbad, California 92008.

### **PRODUCT USE LIMITATION**

This product is developed, designed and sold exclusively *for research purposes and in vitro use only*. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

Please refer [www.thermofisher.com](http://www.thermofisher.com) for Material Safety Data Sheet of the product.

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.