


Pole figure measurement (Schulz reflection method)

[Optic]



[Measurement flow]

A. Test scan

The profile is measured in wide range. Phase identification is performed and sample orientation is investigated compared to intensity ratio of database card. Then, evaluated lattice plane is decided.

B. Investigation of slit condition


Slit condition changes according to sample size and peak width. The peak is measured in lowest alpha angle.

C. Main measurement

Pole figure measurement is performed.

[Exchange of geometry]


1. It is required to finish auto alignment of Focusing method with "Standard sample holder". In most cases, 5deg soller slit is used.
2. Multipurpose attachment is mounted. The power of RCD is turned off. (There is the switch on back side).
3. After connecting attachment cable, the power of RCD is turned on.
4. The sample is mounted on MPA. Rolling direction is faced to beta:0deg.
5. Click "RINT2500 Right system" in [Rigaku]-[Control]-[Rigaku Control Panel]. Select items as follows.



6. Click "Auto alignment" in [Rigaku]-[Control]. The optimization of theta axis and sample height is performed.

[A. Test scan]

- A-1. Click "Standard measurement" in [Rigaku]-[Right measurement].
Enter parameters as follow.




- A-2. Press "Start measurement".

A-3.

Phase identification is performed. And, the intensity ratio of measured peak and database is compared. If there is no difference between measured peak and database card, it seems that there is no orientation in the sample. If there is the difference, the sample is regarded as texture, pole figure measurement is performed.

[B. Investigation of slit condition]

- B-1. Change "Slit" to slit for Pole Figure in Rigaku Control Panel.



- B-2. Set 2mm divergence height limit slit, and mount schulz slit.

- B-3. Click [Pole Figure measurement] in [Rigaku]-[Right Measurement].
Select "None" in [Transmission] and "Schulz reflection" in [Reflection].
Enter sample size, alpha start angle and 2theta angle measured in test scan.

- B-4. Push "Graph of irradiation width" button.
Optimum divergence slit is calculated when pushing "Recalc" button

- B-5. Enter the optimum DS in "DivSlit".

Pole Figure Measurement

File Edit View Option Help

No.1 | No.2 | No.3 | No.4 | No.5 | No.6 | No.7 | No.8 | No.9 | No.10 | No.11 | No.12 | No.13 | No.14 | No.1

Operator: Administrator Memo:

h	k	l	Crystal system	Thickness[cm]	L.absorp.co.[1/cm]	Sample[mm]	<input checked="" type="checkbox"/> Multiple profile
1	1	1	Cubic	0.0300000	470.0000	50.0000000	

Method: Continuous Counting unit: cps Auto Scale Init scale: 1000

Main measurement

Method	Alpha start	Alpha stop	Alpha step	kV	mA	BG meas.
Transmission	None					
Reflection	Schulz reflc	25.000	90.000	5.000	50	300
						1 meas.

Scanning	Beta start	Beta stop	Beta step	Beta scan rate	Gamma width
Coaxial	0.00	360.00	5.0	360.00	0mm

2theta angle	DivH.L.Slit	DivSlit	ScfSlit	RecSlit	Comment
Reflection	38.4100	2mm	1/3deg.	7.00mm	7.00mm

BG.measurement


BG Scanning	Beta start	Beta stop	Beta step	Beta scan rate
Scan.	0.00	360.00	5.0	360.00

Alpha step	BG.1 angle	Scattering S.	Receiving S.	BG.2 angle	Scattering S.	Receiving S.
Transmission						
Reflection	5.000	41.0000	7.00mm	7.00mm		


Pre-measurement

Peak Angle	Peak low ang.	Peak high ang.	BG1 low ang.	BG1 high ang.	BG2 low ang.	BG2 high ang.
Reflection	43.5018	42.5800	44.1600	40.1958	41.8041	


Estimated Measurement Time Transmission:0minutes Reflection:38minutes 11:15




B-6. Click [Manual measurement] in [Rigaku]-[Right Measurement].
 Set [DivSlit] to optimum DS, [RecSlit] to 0.45mm. The value of [SctSlit] is same as [DivSlit].
 Profile measurement is performed.




B-7. Read off peak position.




B-8. Set to the position with "Move".




B-9. Open [Attachment Control] window. Move to the lowest alpha angle to be measured in [Alpha-R]. The max range of Alpha-R is 15deg to 90deg. As alpha angle is lower, the peak width broadens.




B-10. Beta scan is performed to find the strongest point.




B-11. Move to the strongest beta angle.



B-12. 2theta measurement is performed to decide low and high position of peak and BG range.



B-13. Press “peak search” button, left click in low position and right click in high position. The peak position is displayed on the bottom of window.



B-14. Enter “Peak Angle” measured in alpha90deg, low angle and high angle. The software selects the optimum RS automatically from registered RS value in “Rigaku Control Panel”.

There are two types in BG measurement. [Scan] measures specified beta range. However, only beta point in which the intensity is the lowest is saved. [1point] performs BG measurement in beta point that is the lowest intensity in main measurement.

Measurement time depends on the rate of beta scan. In most cases, beta scan rate is 120 to 360 deg/min. Estimated measurement time is displayed on the bottom of window.

Gamma oscillation is used when grain size is large. Although 5, 10 and 20mm can be selected on software, oscillation width depends on the size of cum equipped with MPA in practice

B-15. 1 measurement or 2 measurement can be selected in BG measurement.
 Enter center BG position to be measured in BG.1 or BG.2 angle. Adjust the value not to overlap low and high of peak. In case there is close peak, change RS to smaller. Even if RS value of BG is different from RS of main measurement, software can correct the difference of intensity due to different RS.

Pole Figure Measurement

File Edit View Option Help

No.1 | No.2 | No.3 | No.4 | No.5 | No.6 | No.7 | No.8 | No.9 | No.10 | No.11 | No.12 | No.13 | No.14 | No.15

Operator: Administrator Memo:

h	k	l	Crystal system	Thickness[cm]	L.absorp.co.[1/cm]	Sample[mm]
1	1	1	Cubic	0.0300000	470.0000	50.0000000

Method: Continuous Counting unit: cps Auto Scale Init scale: 1000

Main measurement

Method	Alpha start	Alpha stop	Alpha step	kV	mA	BG meas.	
Transmission	None						
Reflection	Schulz reflec	25.000	90.000	5.000	50	300	2 meas.

Scanning	Beta start	Beta stop	Beta step	Beta scan rate	Gamma width
Coaxial	0.00	360.00	5.0	360.00	0mm

2theta angle	DivH.L.Slit	DivSlit	SctSlit	RecSlit	Comment
Reflection	38.4108	2mm	1/3deg.	7.00mm	7.00mm

BG.measurement


BG Scanning	Beta start	Beta stop	Beta step	Beta scan rate
Scan	0.00	360.00	5.0	360.00

Alpha step	BG.1 angle	Scattering S.	Receiving S.	BG.2 angle	Scattering S.	Receiving S.
Transmission						
Reflection	5.000	36.7500	7.00mm	7.00mm	40.2500	7.00mm

Pre-measurement


Peak Angle	Peak low ang.	Peak high ang.	BG1 low ang.	BG1 high ang.	BG2 low ang.	BG2 high ang.
Reflection	38.4108	37.6800	39.3200	36.0464	37.4536	39.5464

Estimated Measurement Time Transmission:0minutes Reflection:57minutes 12:10



*Defocus correction

It is important to perform “defocus correction” to draw pole figure correctly in reflection method. As alpha angle is smaller, the intensity decreases due to defocus. Non-oriented bulk sample or powder sample are used for defocus correction. Basically, the shape of defocus curve depends on goniometer radius, 2theta angle and the width of RS. DS has little affect on the intensity decrease. Therefore, when measuring non-oriented sample, it need to use the same RS as real sample.



Cubic Crystal Orientation Display program Cubic COD



Miller Indices



h k l
u v w

Euler Angles

ϕ_1 Φ ϕ_2

Position Disp size

Background Color Line size



[operation]

How to open [Cubic COD]---
Doubleclick "CubicCODIsp.jar" in C:\CTR\bin

Miller Indices

h k l
u v w

Euler Angles

$\phi 1$ Φ $\phi 2$

Position

Disp size

Background Color

Line size

(hkl)[uvw] is selected on [Miller Indices].
Euler angle is calculated on [Euler Angles] by pushing "Calc" button.
Program does not calculate Euler angle and all values of Euler angle display "0.0" when $h*u+k*v+l*w$ does not equal to zero.

Euler angle by Bunge type is entered in $\phi 1$, Φ , $\phi 2$.
Calculated Miller Indices is displayed on [Miller Indices] with "Calc" button.
If $h*u+k*v+l*w$ does not equal to zero, all values of Euler angle display "0.0"
All values of Euler angle display "0.0" by "Clear" button.

[Position]:not used.


[Disp size]:size of 3D orientation window.

[Background Color]: exchange of background color in 3D orientation window.

[Line size]: change of line width in 3D orientation window.

"Disp" button: 3D orientation window is displayed. Enable to display multiple windows.

[3D orientation window]



In above window, original point is located at back of left below.
a axis is green, b axis is orange and c axis is magenta.

Purple triangle expresses plane (111)

Direction is expressed with white arrow on purple triangle [1-10]

Movement of lattice : right drag with mouse

Rotate : left drag with mouse

Zoom-in and -out : scroll wheel drag with mouse

[Summary]

It is possible to specify crystal orientation from pole figure to some extent. Or, ODF analysis enables to get Euler angle that indicates the relation between material axis(ND, RD) and lattice plane. It is possible to calculate crystal orientation {hkl}<uvw> from this Euler angle.

However, it is difficult to estimate actual Plane-Direction from {hkl}<uvw>.

This program was made to understand crystal orientation quickly by 3D display.

{hkl}<uvw> means identification direction. (hkl)[uvw] is expressed as singularity expression.

This program uses singularity expression. After this, crystal orientation is expressed as (hkl)[uvw].

[Features]

1. Can calculate Euler angle from (hkl)[uvw].
2. Can calculate (hkl)[uvw] from Euler angle
3. Can display crystal orientation with 3D.
4. Can display (hkl) plane without Euler angle.
Program does not calculate Euler angle when $h*u+k*v+l*w$ does not equal to zero.
In this case, direction is not on plane.
5. It is possible to zoom-in and -out and rotate 3D lattice.

[Angle calculation]

Crystal orientation is expressed as (hkl)[uvw].

Euler angle is expressed as $\phi 1, \Phi, \phi 2$.

$$h=n*\sin\Phi*\sin\phi 2$$

$$k=n*\sin\Phi*\cos\phi 2$$

$$l=n*\cos\Phi$$

$$u=m(\cos\phi 1\cos\phi 2 - \sin\phi 1\sin\phi 2\cos\Phi)$$

$$v= m(-\cos\phi 1\sin\phi 2 - \sin\phi 1\cos\phi 2\cos\Phi)$$

$$w=m*\sin\phi 1\sin\Phi$$

$$\phi 1=\arcsin(w/\sqrt{(u*u+v*v+w*w)}*\sqrt{(\sqrt{(h*h+k*k+l*l))/(h*h+k*k)})}$$

$$\Phi = \arccos(1/\sqrt{(h*h+k*k+l*l)})$$

$$\phi 2=\arccos(k/\sqrt{(h*h+k*k)})=\arcsin(1/\sqrt{(h*h+k*k)})$$

where $h*u+k*v+l*w$ equals to zero.

[operation environment]

It is required that Java runtime and Java 3D are installed.