

Procedure of SEM sample preparation in EM lab. (morphology)

Fix sample in Fixation Buffer* for 4hr in room temp. or overnight in 4°C buffer and rinses 3X 15 min in room temp. before proceed to following steps

Buffer rinses	3X 15 min in room temp.
30% EtOH	10 min
50% EtOH	10 min
70% EtOH	10 min or overnoight
85% EtOH	20 min
95% EtOH	20 min
100% EtOH	20 min
100% EtOH	20 min or overnight
100% Acetone	20 min ^s
100% Acetone	20 min

Critical point drying or freeze drying

Sample on stubs

Coating Au

Observation

Note:

*Fixation Buffers (choose one):

(1) Gluteraldehyde buffer

2.5% GA + 4% PFA/0.1M phosphate or cacodylate buffer

The buffer should be pH to 7.2, store at 4°C (good for one year)

(2) Osmium tetroxide buffer

1% OsO₄/0.1M phosphate or cacodylate buffer 4hr in room temp.

Make fresh - 1-2 days in advance to dissolve

Use FUME HOOD! and use gloves all the time

This buffer can be used after buffer (1)

^sDo not proceed to acetone step if you are not going to the drying procedure right away, as the materials will be fragile if left in acetone for too long.

GA = gluteraldehyde

PFA = paraformaldehyde

OsO₄ = osmium tetroxide