Jer-Ming Hu Feb. 21, 2001

& Microtoming

After the samples are embedded in paraffin, it is ready for microtome dissection. Stick the sample in a wooden block which has been soaked in paraffin before hand. For general dissection of plant tissue, 8-12.5 μ m thick is typical, but varies on different materials. Cut the ribbon in desired length and place on a warm water bath (55°C), this could be a Petri dish or a designed instrument. Use microscope slide to pull out the ribbon and place on the middle of the slide. Dry the slides on a warm plate (55°C) overnight.

& Staining (followed Johansen, 1940)

To remove paraffin, follow the steps and transfer to next step every 5-10 mins.

xylene xylene + ethanol (1:1) 95% EtOH 70% EtOH 35% EtOH dH_2O **Tol-Blue** in dH_2O – or stain in other dye (**Safranin***) dH_2O dH_2O dH_2O 35% EtOH 70% EtOH 95% EtOH ysylene + ethanol (1:1) xylene

Take out the slide and add a drop of Permont, put on the cover glass and put a "weight" on the slide. The slide should be ready overnight.

*Safranin solution (for 200ml)

dissolve 2g **Safranin O** in 100ml **methyl cellosolve** (=ethylene glycol) When solution complete, add 50ml **95%EtOH** and 50ml **dH**₂**O** Then add 2g **NaOAc** and 4ml **formalin**

t <u>Things I use in Mike's lab</u>:

Microscope slide – Fisherbrand superfrost[®]/Plus, precleaned (Cat.#12-550-15) Cover glass – Corning, No. 1, 22mm sq., FISHERfinest[™] Premium cover glass Microtome – Leica RM2065; and Leica disposable microtome blades (Model818) Weight – West Coast premium magnum shot No. 9, West Coast Shot Inc., filled in a small glass vial (Vial, S/T, Type I glass, FISHER Cat.#03-338-25B, \$105.3/144vial), but I think the cheaper one will work too (FISHER Cat.#03-377A \$12.22/100vial) **Toluidine Blue O** – Fisher (Cat.#BP107-10)

Safranin O – Sigma (No. S 8884)

Methyl cellosolve (=ethylene glycol=monoethyl ether) – Fisher (Cat.#E182-500)

Reference:

Johansen, D. A. 1940. Plant microtechnique. McGraw-Hill Book Company. New York.