

***E. coli* Competent cell preparation**

(This protocol is from Erh-Min Lai, based on Molecular Cloning, 2nd ed.)

1. Inoculate fresh *E. coli* colony in 3ml SOB+20mM MgSO₄ broth, incubate at 37°C, 16hr.
2. Transfer 1ml culture into 99 ml SOB+20mM MgSO₄, incubate at 37°C, 2.5-3.0 hr, 測 OD₆₀₀。約 0.35-0.4 即可。
3. Transfer to sterile, ice-cold 50 ml PP tube. Cool the cultures to 0°C by storing the tubes on ice for 10 minutes.
4. Centrifugation at 3000g for 10 minutes at 4°C.
5. Decant the medium from the cell pellets. Stand the tubes in an inverted position for 1 minute to allow the last traces of medium to drain away.
6. Resuspend the pellets by swirling or gentle vortexing in ~20 ml (per 50 ml tubes) of ice-cold TEB (for fresh) or FSB (for frozen) transformation buffer. Store the resuspended cells on ice for 10 minutes.
7. Centrifugation at 3000g for 10 minutes at 4°C.
8. Decant the medium from the cell pellets. Stand the tubes in an inverted position for 1 minute to allow the last traces of medium to drain away.
9. Resuspend the pellets by swirling or gentle vortexing in 4 ml (per 50 ml tubes) of ice-cold TEB (for fresh use) or FSB (for store at -70°C and used at a later date).
10. To prepare fresh competent cells
 - a. Add 140 µl of DnD solution into the center of each cell suspension. Immediately mix by swirling gently, and store the suspension on ice for 15 minutes.
 - b. Add an additional 140 µl of DnD solution to each suspension. Mix by swirling gently, and store the suspension on ice for a further 15 minutes.
 - c. Dispense aliquots of the suspensions into chilled, sterile 17x100 ml PP tubes. Store the tube on ice.
11. To prepare frozen stocks of competent cells
 - a. Add 140 µl of DMSO per 4 ml of suspended cells. Mix gently by swirling, and store the suspension on ice for 15 minutes.
 - b. Add an additional 140 µl of DMSO to each suspension. Mix gently by swirling, and then return the suspensions to an ice bath.
 - c. Dispense 100 µl aliquots of the suspensions into chilled, sterile microfuge tubes or tissue culture vials. 以液態氮迅速冷凍，存放於-70°C使用時再取出。
 - d. 使用時由-70°C取出，以手的溫度融解，融了之後立即放於冰上 10 分鐘。
 - e. 分裝至 17x100 mm 的 PP tube，放置於冰上
12. 每 100 µl competent cells 加 25-100 ng transforming DNA(5µl)(體積勿超過 5%)，每次作 2 個 control，一個加入已知量的 plasmid DNA，另一個不加，置於冰上 30mins。
13. 42°C 水浴 90 秒。
14. 冰浴 1-2 分鐘。
15. 加 800 µl SOC medium，37°C 水浴加熱，移至 37°C 培養 45 分鐘。
16. LB+抗生素的 plate 上加入 100-200 µl, spread evenly(用 sterilized 玻璃珠)。
17. Spread 均勻後倒掉玻璃珠。
18. 37°C 培養 12-16 小時可得 transformed colonies。